

Vacuum Reg. Pirol 2.0 -90,VAC,RAIL,NIST



Manufacturer: Greggersen

Product 904605

Number:

Description

Vacuum Regulator Pirol -90, Rail-Mounted Unit NIST

Medical vacuum regulator for dosing vacuum.

With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel. The adjustment mechanism enables stepless, precise and intuitive regulation of the vacuum. The rotatable vacuum gauge allows the user to read the set value from an optimal viewing angle.

The regulator unit meets the requirements of DIN EN ISO 10079-3 (high vacuum /high flow).

Material:

Housing: Aluminium, anodised

Handwheel: Plastic

Rail clamp: Aluminium, anodised

Technical data:

Gas type: Vacuum / VAC
Inlet: NIST housing according to DIN EN ISO 18082
Control range: 0 to -90 kPa ± 10%
Flow rate: > 20 l Freeflow/min
Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum
----------	------------
