


Hose, AIR, 3m, DIN - DIN Coupling Colour-Neutral



Manufacturer: Greggersen

Product 902535

Number:

Description

Low-Pressure Hose

acc. to DIN EN ISO 5359

Angled Probe with reverse flow restrictor,
hose and coupling acc. to DIN

Specification:

Gas type: Compressed Air / AIR

Inlet: Probe acc. to DIN 13260 Part 2

Hose length: 3.0 m

Colour coding: Colour-Neutral

Outlet: Coupling acc. to DIN 13260 Part 2

Properties

Gas Type AIR Compressed Air