

Kolibri 2 Probe DIN O2 Standard



Manufacturer: Greggersen

Product 906400

Number:

Description

KOLIBRI 2, O₂, 0...15 l/min

- Standard DIN PLUG-IN UNIT -

Flow Settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selection device

Medical compact flowmeter for connection

to a terminal unit, with handwheel - incrementally

adjustable flow rate, outlet 9/16-18 UNF

Setting device ergonomically optimised and readability

of the set flow setting possible from different viewing directions

and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas Type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF thread

Flow Rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	DIN
Variant	Standard