

Compressed-Air Vacuum Regulator Specht -90 NF

Probe



Manufacturer: Greggersen

Product 901512

Number:

Description

Compressed-Air Vacuum Regulator Specht-90, Plug-in Unit NF

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine-dosing and quick-action shut-off valve, silencer and filter holder.

Technical data:

Gas type: Compressed air / AIR

Inlet: Probe according to NF S 90-116

Control range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------