

Vacuum Regulator Spatz -10 Rail NIST



Manufacturer: Greggersen

Product 902675

Number:

Description

Vacuum Regulator Spatz-10, Rail-Mounted Unit NIST

Medical vacuum regulator with vacuum gauge,
fine dosing and quick-action shut-off valve,
for mounting on an equipment rail (25x10 mm).
The gas inlet is integrated into the rail clamp (NIST).

Accessories required for operation:

Low-Pressure Hose VAC with NIST outlet
(not included in scope of delivery)

Technical data:

Gas type: Vacuum / VAC

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -10 kPa (0 to -100 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple