

Kolibri 2 Probe MEDAP O2 Paediatric



Manufacturer: Greggersen

Product 906452

Number:

Description

KOLIBRI 2, O₂, 0...5 l/min

- Paediatric MEDAP PLUG-IN UNIT -

Flow settings:

0 / 0.1 / 0.3 / 0.5 / 0.6 / 0.8 / 1 / 2 / 3 / 4 / 5 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices with selector device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment mechanism and readability
of the set flow setting from various viewing directions
and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: Oxygen / O₂

Connection: MEDAP Probe

Outlet: 9/16-18 UNF

Flow rate: 0...5 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	MEDAP
Variant	Paediatric