

## Kolibri, Probe, DIN, AIR OFF/8



**Manufacturer:** Greggersen

**Product** 901403

**Number:**

### Description

Kolibri, for Compressed Air

-Plug-in Unit-

Medical Compact Flowmeter according to DIN EN ISO 15002,  
for coupling into a Terminal Unit, with handwheel  
incrementally adjustable flow rate,  
Outlet 9/16-18 UNF and Tubing Nipple

Gas Type: Compressed Air

Connection: Probe according to DIN 13260 Part 2

Overall Length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18 UNF/ Tubing Nipple

Flow Rate: 0 / 8 l/min

Flow Settings: 0 + 8

# Properties

|          |                    |
|----------|--------------------|
| Gas Type | AIR Compressed Air |
| Variant  | Air Off/8          |