

Forano, Front Part, DIN, O2, neut



Manufacturer: Greggersen

Product 905110

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral - Version III from 08/2019

Gas type: Oxygen / O2

Inscription: oxygen / O2

Gas-specific socket assembly for mounting
on FORANO base part.

Socket assembly according to DIN 13260 Part 2
with parking and operating position.

Release via release ring.

Colour coding:

Colour-Neutral (black)

Release ring inscription: printed

Language: English

Cartridge incl. sealing elements

Technical data:

Gas type: Oxygen / O₂

Inscription: oxygen / O₂

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed inscription
- cartridge with sealing elements

Properties

Gas Type

O₂-Oxygen
