


# Hose, O2, 5m, DIN - DIN Coupling, Colour-Neutral



**Manufacturer:** Greggersen

**Product** 902530

**Number:**

## Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,  
hose and coupling according to DIN

Specification:

Gas type: Oxygen / O<sub>2</sub>

Inlet: Probe according to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: Coupling according to DIN 13260 Part 2

## Properties

**Gas Type**

O<sub>2</sub>-Oxygen