

Kolibri, Probe, DIN, O2 Standard, extended



Manufacturer: Greggersen

Product 901496

Number:

Description

Kolibri, O2, 0...15 l/min

- PLUG-IN UNIT, extended -

Medical compact flowmeter according to DIN EN ISO 15002,
for connection to a terminal unit, with handwheel
step-adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas type: Oxygen / O2

Connection: Probe according to DIN 13260 Part 2

Total length: Probe front part to centre of flow outlet 133 mm

Outlet: 9/16-18 UNF / Tubing Nipple

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
