

Forano, Base Part, Straight Pipe, DIN, AIR



Manufacturer: Greggersen

Product 903201

Number:

Description

Terminal Unit Forano Base Part (Base Part) DIN/AGA/NF

Gas-specific base part for installation in
supply rails or ducts.

Base part with spacer and mounting bridge;
brazed copper pipe, non-return and maintenance valve integrated.

Scope of delivery:

- gas-specific base part with Cu pipe
- spacer
- mounting bridge

Technical data:

Gas type: Compressed Air / AIR

Installation type: Media supply rail

Inlet: Cu pipe 8x1 mm straight, top

Properties

Gas Type AIR Compressed Air