


Falke with Flow Paediatric O2 DIN 9/16



Manufacturer: Greggersen

Product 910352

Number:

Description

Falke with Flow, Paediatric

Compact pressure regulator with contents gauge;
incrementally adjustable flow rate.

Gas type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (Tubing Nipple) bottom

Flow rate: 0...5 l/min

Flow settings: 0/ 0.1/ 0.3/ 0.5/ 0.6/ 0.8/ 1/ 2/ 3/ 4/ 5 l/min

Properties

Gas Type

O₂-Oxygen