

Forano, Base Part, Straight Pipe, DIN, CO2



Manufacturer: Greggersen

Product 903205

Number:

Description

Terminal Unit Forano Base Block (Base Part) DIN/AGA/NF

Gas-specific base block for installation in supply rails or ducts.

Base block with spacer disc and mounting bracket; brazed copper pipe, integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base block with Cu pipe
- spacer disc
- mounting bracket

Technical data:

Gas type: Carbon dioxide / CO2

Installation type: Media supply rail

Inlet: Cu pipe 8x1 mm straight, top

Properties

Gas Type CO2-Carbon Dioxide