

Vacuum Regulator Spatz -90 DIN Probe



Manufacturer: Greggersen

Product 900900

Number:

Description

Vacuum Regulator Spatz-90, Plug-in Unit DIN

Medical vacuum regulator according to the bypass principle,
with vacuum gauge, fine adjustment and quick-acting shut-off valve.

Technical Data:

Gas type: Vacuum / VAC

Inlet: Probe to DIN 13260 Part 2

Control range: -20 to -85 kPa (-0.2 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple