

Kolibri 2 Probe NF AIR Standard



Manufacturer: Greggersen

Product 906432

Number:

Description

KOLIBRI 2, AIR, 0...15 l/min

- Standard NF PLUG-IN UNIT -

Flow Settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel - step-
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: med. Compressed Air / AIR
Connection: Probe according to NF S 90-116
Outlet: 9/16-18 UNF
Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:
900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	NF
Variant	Standard