


Pressure Regulator Unit, AIR, Size 1, 5 bar



Manufacturer: Greggersen

Product 325810

Number:

Description

Pressure Regulator Unit for reducing vessel pressure to 5 bar, design based on DIN EN ISO 7396-1, each pressure regulator stage duplicated with shut-off valves, safety valve, contact transmitter, fully mounted on a mounting plate.

Size: 1

Flow Rate: 50 Nm³/h

Gas Type: Compressed Air / AIR

Pressure: 500 kPa (5 bar)

Contact Transmitter: 400 / 600 kPa (4 bar / 6 bar)

Properties

Gas Type AIR Compressed Air