

Vacuum Reg. Pirol 2.0 -90,VAC,PROBE,CAR



Manufacturer: Greggersen

Product 904620

Number:

Description

Vacuum Regulator Pirol -90, Plug-in Unit CARBAMED

Medical vacuum regulator for vacuum dosing.

With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel.

The adjustment mechanism enables infinitely variable, accurate and intuitive vacuum regulation. The rotatable vacuum gauge allows the user to read the set value from an optimum viewing angle.

The PIROL regulator units meet the requirements of EN ISO 10079-3 (high vacuum / high flow).

Material:

Housing: Aluminium, anodised

Handwheel: PA plastic

Probe: Brass, chrome-plated

Technical data:

Gas type: Vacuum / VAC

Inlet: CARBAMED probe according to SN ENV 737-6:2005

Control range: 0 to -90 kPa $\pm$ 10%

Flow rate: > 20 l freeflow/min at -90 kPa

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum
-----------------	------------
