

Vacuum Reg. Skua 2.0 -90,AIR,RAIL,NIST



Manufacturer: Greggersen

Product 904601

Number:

Description

Compressed-Air Vacuum Regulator Skua -90, Rail-Mounted Unit NIST

Medical compressed-air vacuum regulator for dosing vacuum.

The vacuum is generated by means of compressed air according to the Venturi principle. With 360° rotatable vacuum gauge; large, ergonomic, adjustment handwheel; rotatable exhaust sleeve.

The gas inlet is integrated into the rail clamp (NIST).

The adjustment mechanism allows for infinitely variable, precise and intuitive vacuum regulation. The rotatable vacuum gauge enables the user to read the set value from an optimal viewing angle.

The air flow can also be aligned individually via the exhaust sleeve.

The regulator unit complies with the requirements of DIN EN ISO 10079-3:2014 (high vacuum /high flow).

Material:

Housing: aluminium, anodised

Handwheel: plastic

Rail clamp: aluminium, anodised

Technical data:

Gas type: med. compressed air / AIR

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -85 kPa $\pm$ 10%

Flow rate: >20 l Freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type

AIR Compressed Air
