

Hose, CO2, 1.5m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900045

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with backflow limiter,
hose and equipment connection NIST

Specification:

Gas type: Carbon dioxide / CO2

Inlet: Probe according to DIN 13260 Part 2

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type CO2-Carbon Dioxide