


Falke w. Flow Neo O2 DIN Tubing Nipple



Manufacturer: Greggersen

Product 910203

Number:

Description

Falke with Flow, NEO,
Compact pressure regulator with contents gauge;
incrementally adjustable flow rate.

Gas type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing Nipple bottom

Flow rate: 0...1 l/min

Flow Settings: 0/ 0.1/ 0.15/ 0.2/ 0.25/ 0.3/ 0.4/ 0.5/ 0.6/ 0.8/ 1.0 l/min

Properties

Gas Type O₂-Oxygen