


Angled Probe VAC Hose Barb 7mm, DIN, ISO V4A



Manufacturer: Greggersen

Product 505152

Number:

Description

Angled Probe for Vacuum

with 7 mm hose barb

Probe front part made of 1.4404 (stainless steel)

Gas type: Vacuum / VAC

Connection: Probe according to DIN 13260 Part 2

Outlet: 7 mm Hose Barb

Colour coding: according to ISO 9170-1 (yellow)

Properties

Gas Type	VAC-Vacuum
-----------------	------------