

Vacuum Regulator Spatz -90 Probe DIN ANTIM



Manufacturer: Greggersen

Product 500243

Number:

Description

Vacuum Regulator Spatz-90, Plug-in Unit DIN, non-magnetic

Medical vacuum regulator based on the bypass principle,
with vacuum gauge, fine-adjustment and quick-shut-off valve.

Technical data:

Gas type: Vacuum / VAC

Inlet: Probe to DIN 13260 Part 2

Control range: -20 to -85 kPa (-0.2 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple