


Forano, Front Part, AGA, AGSS (UTL), ISO



Manufacturer: Greggersen

Product 905147

Number:

Description

Terminal Unit FORANO - Front Part AGA

Gas type: AGSS (UTL)

Labelling: AGSS / UTL

Gas-specific Socket Assembly for mounting
on Base Part FORANO. Design complies with
DIN EN ISO 9170-1;

Socket Assembly according to SS 87 524 30

Release via Release Ring.

Colour coding according to ISO 32.

Release Ring labelling: printed

Cartridge incl. sealing elements

Technical data:

Gas type: AGSS (UTL)

Labelling: AGSS / UTL

Colour coding: ISO 32

Front diameter: Ø 42mm

Scope of delivery:

- Gas-specific Socket Assembly
- Release Ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type	AGSS-UTL
-----------------	----------
