


Falke without Flow N2O DIN +1DIN (small cylinder)



Manufacturer: Greggersen

Product 910379

Number:

Description

Falke without Flow with a DIN quick coupling,
Compact Pressure Regulator with contents gauge,
DIN quick coupling at a pressure of 450 kPa.

Gas type: Nitrous oxide / N2O

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 12 (small cylinder)

Outlet: Quick coupling DIN 13260 Part 2, top

Output: 450 kPa

Properties

Gas Type N2O-Nitrous Oxide