


Hose, O2, 3m, DIN - DIN Coupling, cn



Manufacturer: Greggersen

Product 902534

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled probe with reverse flow limiter,
hose and coupling according to DIN

Specification:

Gas type: Oxygen / O₂

Inlet: Probe according to DIN 13260 Part 2

Hose length: 3.0 m

Colour coding: Colour-Neutral

Outlet: Coupling according to DIN 13260 Part 2

Properties

Gas Type

O₂-Oxygen