


Kolibri, Rail AIR Standard - NIST



Manufacturer: Greggersen

Product 902608

Number:

Description

Kolibri, AIR, 0...15 l/min

RAIL-MOUNTED UNIT NIST

Medical Compact Flowmeter

according to DIN EN ISO 15002; flow is set via

handwheel in 10 flow settings;

gas inlet NIST integrated into rail clamp;

outlet 9/16-18 UNF and tubing nipple

Gas Type: Compressed Air AIR

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF / Tubing Nipple

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	AIR Compressed Air
----------	--------------------
