

Kolibri 2 Probe DIN AIR Standard+ long



Manufacturer: Greggersen

Product 906412

Number:

Description

KOLIBRI 2, AIR, 0...30 l/min

- Standard+ DIN PLUG-IN UNIT -

Long version

Flow settings:

0 / 2 / 3 / 4 / 6 / 8 / 12 / 15 / 20 / 25 / 30 l/min

KOLIBRI 2 Compact Flowmeter

with fixed orifices and selector mechanism

Overall length "long"

Probe front part to centre of flow outlet

(Distance from terminal unit to tubing nipple):

75 mm (intermediate piece 60 mm)

Medical compact flowmeter for connection

to a terminal unit, with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> increased safety when
setting the correct flow rate

Gas type: Med. Compressed Air / AIR
Connection: Probe according to DIN 13260 Part 2
Outlet: 9/16-18 UNF
Flow rate: 0...30 l/min at inlet pressure of 450 kPa (4.5 bar)

Accessories required for operation:
900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	DIN
Variant	Standard+ Long