


Forano, Front Part, DIN, O2 93, ISO



Manufacturer: Greggersen

Product 905108

Number:

Description

Terminal Unit FORANO - Front Part DIN

ISO 32 -

Gas type: Oxygen 93 / O2 93

Labelling: oxygen 93 / O2 93

Gas-specific socket assembly for mounting
on base part FORANO.

Socket assembly according to DIN 13260 Part 2
with park and operating position. Release via release ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - white)

Release ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas type: Oxygen 93 / O2 93

Labelling: oxygen 93 / O2 93

Colour coding: according to ISO 9170-1 (white)

Front diameter: Ø 42mm

Scope of Delivery:

- Gas-specific socket assembly
- Release ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type	O293-Oxygen 93
-----------------	----------------
