

Hose, O2, 1.5m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900038

Number:

Description

Low-Pressure Hose

in accordance with DIN EN ISO 5359

Angled Probe with reverse flow restrictor,
hose and equipment connection NIST

Specification:

Gas type: Oxygen / O2

Inlet: Probe in accordance with DIN 13260 Part 2

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type

O2-Oxygen