

Hose, O2, 0.55m, DIN - NIST, fn



Manufacturer: Greggersen

**Product
number:** 500390

Description

Low-Pressure Hose
according to DIN EN ISO 5359

Angled probe with non-return valve,
hose and NIST equipment connection

Specification:

Gas type: Oxygen / O2

Inlet: Probe according to DIN 13260 Part 2

Hose length: 0.55 m

Colour code: Colour-neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type O2-Oxygen