

Kolibri 2 Probe MEDAP AIR Standard



Manufacturer: Greggersen

Product 906454

Number:

Description

KOLIBRI 2, AIR, 0...15 l/min

- Standard MEDAP PLUG-IN UNIT -

Flow Settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Adjustment device ergonomically optimised and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: Med. Compressed Air / AIR

Connection: MEDAP Probe

Outlet: 9/16-18 UNF

Flow Rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose Connection 9/16" Hose Barb - Compressed Gas

Properties

Gas Type	AIR Compressed Air
Inlet Standard	MEDAP
Variant	Standard