

Falke with Flow Standard O2 AGA Tubing Nipple



Manufacturer: Greggersen

Product 910328

Number:

Description

Falke with Flow, Standard,
Compact pressure regulator with high-pressure cylinder connection
to SS 367615, contents gauge, incrementally adjustable
flow rate.

Gas type: Oxygen / O2

Connection: High-pressure cylinder connection to SS 367615; No. E1 (W21.8)

Outlet: Tubing nipple bottom

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O2-Oxygen