

Vacuum Regulator Spatz -90 Rail NIST



Manufacturer: Greggersen

Product 902622

Number:

Description

Vacuum Regulator Spatz-90, Rail-Mounted Unit NIST

Medical vacuum regulator based on the bypass principle,
with vacuum gauge, fine adjustment and quick-shut-off valve,
for attachment to an equipment rail (25x10 mm).
The gas inlet is integrated into the rail clamp (NIST).

Accessories required for operation:

Low-pressure hose VAC with NIST outlet
(not included in delivery)

Technical data:

Gas type: Vacuum / VAC

Inlet: NIST housing according to DIN EN ISO 18082

Control range: -20 to -85 kPa (-0.2 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple