

Vacuum Regulator Spatz -10 DIN Probe



Manufacturer: Greggersen

Product 901915

Number:

Description

Vacuum Regulator Spatz-10, Plug-in Unit DIN

Medical vacuum regulator with large vacuum gauge (Ø 63mm), fine dosing and quick-action shut-off valve.

Technical Data:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -10 kPa (0 to -100 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple