


Forano, Front Part, BS, VAC, ISO



Manufacturer: Greggersen

Product 905162

Number:

Description

Terminal Unit FORANO - Front Part BS

Gas Type: Vacuum / VAC

Labelling: vacuum / VAC

Gas-specific socket assembly for mounting
on FORANO base part.

Socket assembly according to BS 5682.

Release via release ring.

Colour coding:

Socket assembly coloured according to DIN EN ISO 9170-1 (ISO32).

Release ring labelling: printed

Cartridge incl. sealing elements

Technical data:

Gas type: Vacuum / VAC

Labelling: vacuum / VAC

Colour coding: according to DIN EN ISO 9170-1 (ISO32)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed labelling
- cartridge with sealing elements

Properties

Gas Type	VAC-Vacuum
-----------------	------------
