

Nebuliser Plug-In Valve with Socket AIR/O2



Manufacturer: Greggersen

Product 900765

Number:

Description

Fitis Nebuliser Plug-in Valve with Socket AIR/O2

Plug-in Unit

fixed setting with socket (top) and

Combi-Probe, for coupling into a Terminal Unit.

Technical Data:

Gas type: Compressed Air / Oxygen

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow Rate: 5 l/min $\pm$ 20%