

Falke w. Flow Stand. O2 DIN Tubing Nipple +2DIN



Manufacturer: Greggersen

Product 910208

Number:

Description

Falke with Flow, Standard
with 2 quick-release couplings,
compact pressure regulator with contents gauge;
incrementally adjustable flow rate,
DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing nipple / 2 quick-release couplings according to DIN 13260 Part 2

Flow rate: 0...15 l/min tubing nipple - quick-release coupling max. 120 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O₂-Oxygen