

Hose, O2, 1.5m, MEDAP - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 500153

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe,

Hose and Device Connection NIST

Specification:

Gas Type: Oxygen / O2

Inlet: MEDAP

Hose Length: 1.5 m

Colour Coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type

O2-Oxygen