

Vacuum Regulator High-Spatz -90 Rail NIST



Manufacturer: Greggersen

Product 902623

Number:

Description

Vacuum Regulator High-SPATZ-90, VAC, Rail NIST

- RAIL-MOUNTED UNIT -

Gas type: Vacuum

Output: adjustable between 0 and -85 kPa (0 to -0.85 bar)

Medical vacuum regulator, with vacuum gauge,
fine metering and quick-action shut-off valve,
for mounting on an equipment rail (25x10mm)

Accessories required for operation:

Low-Pressure Hose VAC with NIST outlet
(not included in scope of delivery)

RAIL-MOUNTED UNIT

Gas type: Vacuum

Connection: NIST housing according to DIN EN ISO 18082

Outlet: Vacuum Tubing Nipple

Output: adjustable between 0 and -85 kPa (0 to -0.85 bar)