

Forano, Front Part, DIN, N2, Colour-Neutral



Manufacturer: Greggersen

Product 905106

Number:

Description

Terminal Unit FORANO - Front Part DIN

Neutral - Version III from 08/2019

Gas type: Nitrogen / N2

Inscription: nitrogen / N2

Gas-specific Socket Assembly for mounting
on Base Part FORANO.

Socket Assembly according to DIN 13260 Part 2
with park and operating position. Release via Release Ring.

Colour coding:

Colour-Neutral (black)

Release Ring inscription: printed

Language: English;

Cartridge incl. sealing elements;

Technical data:

Gas type: Nitrogen / N2

Inscription: nitrogen / N2

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of delivery:

- Gas-specific Socket Assembly
- Release Ring incl. printed inscription
- Cartridge with sealing elements

Properties

Gas Type	N2-Nitrogen
----------	-------------
