


Hose, AIR, 2m DIN-NIST, Colour-Neutral angled



Manufacturer: Greggersen

Product 505097

Number:

Description

Special Low-Pressure Hose

++ MukoStar ++

Hose for supplying the device

MukoStar with compressed air 5 bar from a terminal unit;

Angled Probe with hose

with Non-Return Valve, hose diameter 6.7mm

and device connection NIST - 90° angled

Specification:

Gas type: Medical Air / AIR

Inlet: Probe according to DIN 13260 Part 2

Hose length: 2 m

Colour coding: Colour-Neutral

Outlet: gas-specific screw connection NIST

90° angled

Properties

Gas Type	AIR Compressed Air
----------	--------------------
