

Compressed-Air Vacuum Regulator Specht -16 Rail NIST



Manufacturer: Greggersen

Product 902626

Number:

Description

Compressed-Air Vacuum Regulator Specht-16, Rail-Mounted Unit NIST

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine dosing and quick-action shut-off valve, silencer and filter holder.

For mounting on an equipment rail (25x10 mm).

The gas inlet is integrated into the rail clamp (NIST).

Accessories required for operation:

Low-pressure hose AIR with NIST outlet

(not included in the scope of delivery)

Technical data:

Gas type: Compressed air / AIR

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -16 kPa (0 to -160 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------
