

Hose, AIR/O₂, 5m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900472

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow restrictor,
hose and NIST device connection

Specification:

Gas type: Compressed Air/Oxygen / AIR/O₂

Inlet: Probe according to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection