


## Forano, Base Part, DVE, DIN, AIR



**Manufacturer:** Greggersen

**Product** 903231

**Number:**

### Description

Terminal Unit Forano Base Part (Base Part) DIN/AGA/NF

Gas-specific base part with hose barb

for installation in ceiling-mounted supply systems.

Base part with mounting bracket, integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with hose barb
- mounting bracket

Technical data:

Gas type: Compressed air / AIR

Installation type: DVE

Inlet: Hose barb 7 mm, rear

### Properties

Gas Type

AIR Compressed Air

---