

Vacuum Reg. Pirol 2.0 -90,VAC,PROBE,DIN



Manufacturer: Greggersen

Product 904604

Number:

Description

Vacuum Regulator Pirol -90, Plug-in Unit DIN

Medical vacuum regulator for vacuum dosing.

With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel. The adjustment mechanism enables stepless, precise and intuitive vacuum regulation. The rotatable vacuum gauge allows the user to read the set value from an optimal viewing angle.

The regulator unit meets the requirements of DIN EN ISO 10079-3 (high vacuum / high flow).

Material:

Housing: aluminium, anodised

Handwheel: plastic

Probe Front Part: stainless steel

Technical Data:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -90 kPa $\pm$ 10%

Flow rate: > 20 l free flow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum
----------	------------
