

Nebuliser Plug-in Valve with Clamp AIR



Manufacturer: Greggersen

Product 500370

Number:

Description

Fitis Nebuliser Plug-in Valve AIR

Plug-in Unit

Preset nebuliser plug-in valve, with side clamp,
for coupling into a terminal unit,
with Tubing Nipple

Technical Data:

Gas type: Compressed Air / AIR

Connection: Probe to DIN 13260 Part 2

Outlet: 9/16-18 UNF/Tubing Nipple

Flow Rate: 5 l/min $\pm$ 20%

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------