

Hose, O2, 3m, DIN - NIST, ISO



Manufacturer: Greggersen

Product 902411

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and device connection NIST

Specification:

Gas type: Oxygen / O₂

Inlet: Probe according to DIN 13260 Part 2

Hose length: 3 m

Colour coding: ISO 32 (white)

Outlet: gas-specific screw connection NIST

Properties

Gas Type

O2-Oxygen
