


Falke without Flow N2O DIN +1DIN (large cylinder)



Manufacturer: Greggersen

Product 910385

Number:

Description

Falke, without flow with 1 quick coupling,
Compact pressure regulator with contents gauge;
with one quick coupling at a pressure of 450 kPa.

Gas type: Nitrous oxide / N2O

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 11 (large cylinder)

Outlet: Quick coupling according to DIN 13260 Part 2, bottom

Output: 450 kPa

Properties

Gas Type N2O-Nitrous Oxide