

Hose, VAC, 3m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900361

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe,

Hose and Equipment Connection NIST

Specification:

Gas Type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Hose Length: 3.0 m

Colour Coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type

VAC-Vacuum