

Pressure Regulator Habicht with Kolibri



Manufacturer: Greggersen

Product 500046

Number:

Description

O2 Pressure Regulator Unit HABICHT

Suction and insufflation unit,

with cylinder connection

Combi unit with cylinder connection, pressure regulator

with high-pressure cylinder connection according to DIN 477, preset

to 500 kPa (5 bar), with contents pressure gauge, safety valve,

ejector with vacuum gauge; Kolibri Standard for

dosing the flow rate;

Scope of delivery without overflow or humidifier bottle

Gas type: Oxygen / O2

Connection: High-pressure cylinder connection according to DIN 477, no. 9

Outlet: Tubing nipple / vacuum nipple

Output: 0...15 l/min; adjustable between 0 and -85 kPa (0 to -0.85 bar)

SN 500046 to SN 500046